

## Tracoe Twist Plus – REF 521-X

### Tracoe Twist Plus spare inner cannulas

#### Product Description

The Tracoe Twist Plus tracheostomy tubes are used as ventilation and therapy tubes with inner cannulas in hospitals, rehabilitation, and community settings. A special feature is the curved neck flange that moves in both the vertical and horizontal axes, following the movements of the patient's head. The thin tube walls of inner and outer cannulas allow for a good inside - outside diameter ratio. Tracoe Twist and Twist Plus tracheostomy tubes are made of polyurethane (PU). The low-pressure cuff is made of medical grade PVC.



#### Scope of Supply

- 10 inner cannulas with 15 mm connector

#### General Information

Product line	Twist Plus
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	n. A.
Handling	
Shelf life	5 years, from date of production
Maximum Useful Life	5 weeks, from first use

#### Packaging

Packaging unit	10 units
Weight	See Table
Dimensions (W x H x L)	125 x 70 x 271
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

#### Material Specifications

Material Outer cannula	Polyurethane
Material Neck Flange	n. A.
Material Cuff	n. A.
Biocompatibility	EN ISO 10993-1

#### Storage



#### Certification

Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K961449

## Product Range

Size	Order Number	ID-IC	BA	Gross Weight	Global Trade Item Number
		mm	°	g	
07	REF 521-07-X	7,0	100	93	4035324018239
08	REF 521-08-X	8,0	100	94	4035324018246
09	REF 521-09-X	9,0	100	98	4035324018253
10	REF 521-10-X	10,0	100	98	4035324018260

## Legend

- ID-IC Inside diameter at bottom of inner cannula
- BA Bending angle



Tracoe Medical GmbH  
Reichelsheimer Str. 1 / 3  
55268 Nieder-Olm  
Germany

Atos Medical AB  
Hyllie Boulevard 17  
215 32 Malmö  
Sweden  
[www.atosmedical.com](http://www.atosmedical.com)