

Tracoe Occlusion Cap– REF 516

Occlusion cap for attachment to any tracheostomy tube with 15 mm connector

Product Description

Speech is our most important means of communication, so losing one's voice after going through a tracheostomy highly affects the patient's life. The Tracoe Essentials products offer therapeutic options, e.g. for weaning from ventilation, stimulation of the oropharyngeal region, and to produce voice.



General Information		Material Specifications	
Product line	Tube accessories	Material	Polyethylene
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm	Biocompatibility	EN ISO 10993-1
Country of origin	Germany		
Handling		Storage	
Shelf life	5 years, as from date of production		
Packaging		Certification	
Packaging unit	5 pieces	Product classification acc. MDR	Is
Weight	18 g	CE Marking	
Dimensions (W x H x L)	144 x 40 x 80 mm	GMDN	47119
Material outer packaging:	Cardboard box	Quality system certified to	EN ISO 13485
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil	FDA registration	Clearance Number K961449 / E285787