

Tracoe Speaking Valves – REF 650-TO-C

Tracoe Phon Assist I Orange - Speaking valve in orange signal color for increased distinguishability from other attachments

Product Description

Speech is our most important means of communication, so losing one's voice after going through a tracheostomy highly affects the patient's life. Tracoe Phon Assist I speaking valve has infinitely adjustable airway resistance and offers therapeutic options, e.g. for weaning from ventilation, stimulation of the oropharyngeal region, and to produce voice.



General Information	
Product line	Speaking Valves
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm
Country of origin	Germany
Handling	
Shelf life	5 years, as from date of production
Maximum useful life	6 months

Packaging	
Packaging unit	5 pieces
Weight	41 g
Dimensions (W x H x L)	272 x 49 x 113 mm
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications	
Material Flap	Silicone
Material housing	Polypropylene
Biocompatibility	EN ISO 10993-1

Storage	

Certification	
Product classification acc. MDR	Is
CE Marking	
GMDN	36071
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K043160 / D020866

Dimensions	
Height	38mm
Diameter (machine end/ connector patient end)	26mm/ 15mm
Length (oxygen supply port)	10 mm



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com